


## 8'-6" CONTROL ROOM **EQUIPMENT** ROOM CUSTOMER TO PROVIDE CLEAR ACCESS FOR SERVICING OF CHILLER REMOTELY LOCATED BY CUSTOMER 8'-8" \_\_\_\_\_\_BARRIER TBD\_\_\_\_\_\_\_\_ PROPOSED MRI EQUIPMENT LAYOUT x1 (THIS INSTALLATION WILL REQUIRE MAGNETIC SHIELDING TO CONTAIN THE 5 GAUSS FIELD WITHIN PROCEDURE ROOM WALLS) SCALE





(THIS INSTALLATION WILL REQUIRE MAGNETIC SHIELDING TO CONTAIN THE 5 GAUSS FIELD WITHIN PROCEDURE ROOM WALLS)

# 

OPART OPEN MRI SYSTEM

#### **GENERAL NOTES**

ALL SPECIFICATIONS AND SPACE REQUIREMENTS ARE SHOWN GENERICALLY. CONTACT GAMMA TECH, INC. FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL EQUIPMENT COMPONENTS AND ADDITIONAL ROOM PREPARATION SPECIFICATIONS.

RECOMMENDED FINISHED CEILING HEIGHT IS 8'-6". AMBIENT ROOM TEMPERATURE: 68° TO 74° F.

HUMIDITY: 40% TO 60% (NON-CONDENSING). ALLOWABLE TEMPERATURE CHANGE: < 5.4° F OVER AN 8 HOUR PERIOD A MINIMUM OF 10 AIR CHANGES PER HOUR IS SUGGESTED

TYPICAL 12KVA STEP-UP TRANSFORMER

KLEEN-RITE, ARCTIC-CHILL OR EQUIVALENT

DIESEL GENERATOR (ONAN OR EQUIVALENT)

FOR REFRIGERATOR

TYPICAL 38280 EXTERNAL

1400 DIESEL GENERATOR

TRANSFER SWITCH

REMOTE INDICATOR PANEL

MAIN INCOMING POWER SUPPLY: 480VAC, 60AMP OR 208VAC, 125AMP, THREE PHASE, 30kVA, DELTA OR WYE CONNECTED. © PLANNING DIVISION

FILE NAME:

TOSHIBA OPART

**REVISION DATE:** 05.14.2009

19.00 9.75

26.00 39.90

12.12 27.00

12.00 4.00 14.00

CONTROL ROOM EQUIPMENT MUST REMAIN OUTSIDE 2 GAUSS FIELD AND EQUIPMENT ROOM MUST REMAIN OUTSIDE 5 GAUSS FIELD AFTER ROOM IS COMPLIANCE WILL REQUIRE MAGNETIC SHIELDING.

### **EQUIPMENT SCHEDULE**

#### EOLIDMENT IDENTIFICATION NUMBER

EQUIPMENT IDENTIFICATION NUMBER						
	BTU / HOUR		,	WEIGHT	(LBS.)	
		DESCRIPTION	DIN w	/ENSION	S (IN.)	
1	5474	TOSHIBA OPART POTENZA-PARAGON MRI .35 TESLA MAGNET	70.00	83.00	72.00	27150
2	615	TOSHIBA OPART MRI MOTORIZED PATIENT COUCH	27.50	179.00	32.00	1092
3	512	TOSHIBA OPART MRI OPERATOR CONSOLE (MONITOR & KEYBOARD)	19.00	18.00	19.00	95
4	12630	TOSHIBA OPART MRI OIE CONTROL CABINET	44.25	35.00	67.75	1043
5	1707	TOSHIBA OPART MRI REFRIGERATOR	20.00	22.00	28.00	220
6	-	TOSHIBA OPART MRI EMERGENCY RAMP DOWN UNIT				
7	682	TOSHIBA OPART MRI INTERFACE PANEL	35.00		57.00	
8	-	TOSHIBA OPART MRI INNERVISION (17" MONITOR, KEYBOARD & PC)	28.00	16.90	26.60	98
OPTIONAL EQUIPMENT LISTED BELOW IS DONE AS A COURTESY. SITE SPECIFIC DEVICES TO BE DETERMINED BY END USER. ENGINEERS OR CONTRACTORS MAY FURNISH AND INSTALL DIFFERENT SITE SPECIFIC DEVICES AT THEIR DISCRETION. EXACT LOCATION TO BE DETERMINED BY ENGINEER OF RECORD. CONSULT OEM SPECIFICATIONS AND REQUIREMENTS.						
9	34500	KLEEN-RITE, ARCTIC-CHILL OR EQUIVALENT EXT. VERTICAL WATER CHILLER	32.00	55.00	60.00	1010
10	5389	TYPICAL VOLTAGE REGULATOR DISTRIBUTION UNIT	30.13	20.75	58.00	850