



TYPICAL ROOM AND EQUIPMENT LAYOUT

GENERAL NOTES

- ALL SPECIFICATIONS AND SPACE REQUIREMENTS
 ARE SHOWN GENERICALLY. CONTACT GAMMA TECHNOLOGIES
 FOR SPECIFIC REQUIREMENTS PERTAINING TO OPTIONAL
 EQUIPMENT COMPONENTS AND ADDITIONAL ROOM
 PREPARATION SPECIFICATIONS.
- LINE VOLTAGE: 208 240V~, 50 OR 60Hz SINGLE PHASE
- ROTATION MOTORIZED, RANGE 135°...+180° VERTICAL MOVEMENT MOTORIZED RANGE - 29.53"...53.15" SID: 25.59" (650MM)

EQUIPMENT SCHEDULE

EQUIPMENT IDENTIFICATION NUMBER

| | BTU / HOUR | | WEIGHT (LBS.) | | | |
|-----|------------|---|---------------|-----------|-------|-----|
| | | DESCRIPTION | | ENSIONS (| ` ' | |
| (1) | 853* | PLANMED SOPHIE MAMMOGRAPHY GANTRY (W/ MAXVIEW POSITIONING SYSTEM) | 59.00 | 35.50 | 84.30 | 352 |
| (a) | | PLANMED RADIATION PROTECTION SCREEN (OPTIONAL) | 30.50 | 34.10 | 74.00 | 132 |

GAMMA TECH, INC



129 Wild Berry Lane Hampstead, NC 28443 Ph (910) 270-9787 Fax (910) 270-9767 www.gammatech.net © PLANNING DIVISION

FILE NAME: PLANMED SOPHIE MAMMOGRAPHY SYSTEM

REVISION DATE: 02/12/03

* DURING STANDBY



SOPHIE WITH MAXVIEW
BREAST POSITIONING SYSTEM